DUROLANE® INSTRUCTIONS FOR USE

Each mL contains:
Hydronic acid stabilized 20 mg
P.S. sodium chloride solution, pH 7 q.s.

Description
DUROLANE is intended to be used for intra-articular injection in mild to moderate traumatic or degenerative hip osteoarthritis. Additionally, DUROLANE is intended to be used for intra-articular injection for symptomatic treatment of mild to moderate osteoarthritis of isolated synovial joints, and for pain following arthroscopic procedures in the presence of mild to moderate osteoarthritis. It should be injected by an authorized physician, or in accordance with local legislation.

DUROLANE contains 20 mg/mL of stabilized non-animal hyaluronic acid in buffered physiological sodium chloride solution pH 7. DUROLANE is a sterile, transparent viscous gel supplied in a 3 mL glass syringe. The product is free of preservatives.

Hyaluronic acid is identical in all living organisms. It is a natural polysaccharide that is present throughout the tissues of the body, with particularly high concentrations in the synovial fluid and the cartilage. DUROLANE is composed of biosynthetically produced hyaluronic acid which has been purified and stabilized. DUROLANE follows the natural metabolic pathway as endogenous hyaluronic acid.

Mode of Action
The body's hyaluronic acid constitutes a natural part of the articular cartilage and is essential for the normal function of cartilage and ligaments and as a shock absorber. Injections of hyaluronic acid in the joints of patients with osteoarthritis aim to diminish the pain and improve the mobility of the joint.

Dosage
DUROLANE is a single injection, single dose preparation and should only be injected once per treatment course. The recommended dose is 3 mL per joint, hip or shoulder joint. The recommended dose is 1-2 mL for small synovial joints (e.g. thumb).

Indications
Symptomatic treatment of mild to moderate hip osteoarthritis. In addition, DUROLANE is also indicated for the symptomatic treatment of osteoarthritis of the knee or moderate osteoarthritis in the ankle, shoulder, elbow, wrist, fingers, thumb and hand. DUROLANE is also indicated for the symptomatic treatment of osteoarthritis within three months of the procedure.

Contraindications
None.

Warnings
DUROLANE should not be injected if the synovial joint is infected or severely inflamed.

DUROLANE should not be injected if there is an active skin infection or sensitive to any of the ingredients of DUROLANE. DUROLANE should not be injected intracapsularly or extra-articularly or in the synovial tissues or capsule.

Do not reutilize DUROLANE as this may damage the syringe.

Precautions
DUROLANE should be used with caution in patients with venous or lymphatic stasis present in the leg.

DUROLANE has not been tested in pregnant or lactating women. It is not indicated for use in children.

A separate syringe of DUROLANE must be used for each individual joint to be treated.

As with any invasive joint procedure there is a small risk of infection.

DUROLANE should not be injected if the patient is known to be sensitive to hyaluronic acid based products.

Local anaesthetics should not be used if the patient is known to be allergic or sensitive to local anaesthetics.

Injection under fluoroscopic control and with the use of a contrast medium should not be made if the patient is known to be allergic or sensitive to the contrast medium.

In clinical studies, reactions in the knee have not been related with a short interval between first and second injection within 6 months.

Increase in injection pressure may indicate incorrect extra-articular placement of the needle. If this occurs, the injection should be stopped immediately and the needle removed.

The effectiveness of DUROLANE following arthroscopic procedures for diagnosis or examination purposes only or in absence of confirmed osteoarthritis of the joint has not been studied.

DUROLANE should be used with caution in patients with pre-existing chondrocalcoses as injection may lead to an increase in the size of the condition.

Adverse Events
The majority of the reported adverse reactions in clinical studies of the knee and hip were described as transient swelling and/or stiffness of the joint. In the majority of cases, adverse reactions were of mild or moderate intensity and only occasionally required treatment with paracetamol or NSAIDs.

The use of other hyaluronic acid preparations in other joints has not involved any additional safety concerns.

None of the other adverse reactions that have been reported are interpreted as acute inflammatory arthritis or allergic reactions and they did not require medical attention in the form of physical intervention, systemic or intracapsular antibiotics.

Interactions
The safety and effectiveness of DUROLANE concomitantly with other intra-articular injectables have not been established.

Additional administration information
DUROLANE should only be injected by an authorized physician (or in accordance with local legislation), familiar with the technique for the intra-articular injection intended to be treated, and in facilities well suited for intra-articular injections.

DUROLANE should be injected using strict aseptic technique.

DUROLANE should be injected into the joint cavity only.

Intra-articular injection in certain synovial joints will require image guidance to ensure accurate placement and avoid extra-articular placement of the injection.

The route for intra-articular injection with or without image guidance should be chosen so that damage to adjacent tissues is avoided.

The injection site should be swabbed with alcohol or other suitable antiseptic solution before injection.

Remove joint effusion, if present, before injecting for both knee and shoulder. This should be done for both removal of effusion and injection of DUROLANE.

The recommended needle size is 18 to 22 G and with adequate length.

Image guidance becomes increasingly required to deliver the product.

Additional information for treatment of synovial joints requiring image guidance
The intra-articular injection in the hip joints should be given image guidance to ensure accurate placement and avoid extra-articular placement of the injection.

Guidance of other synovial joints is at the discretion of the treating physician.

Injection discomfort can be minimized by use of topical freezing agents or subcutaneously delivered local anaesthetics.

Image guided injection should only be performed if physicians experienced in this type of administration.

Additional information for treatment post arthroscopy
Following the arthroscopic procedure, intra-articular injection should be performed outside the sterile field as the exterior syringe is not sterile.

Joint pain typically under arthroscopic procedures are the knee, hip, shoulder, elbow, ankle, and wrist joints.

Please inform your patient that:
As with any invasive joint procedure it is recommended that patients avoid strenuous activities (e.g. tennis, jogging or rowing) for the first two days after the injection.

Some transient reactions related to the injection of DUROLANE, such as pain and/or redness at the injection site, are to be anticipated. If the symptoms last for more than a week a physician should be contacted.

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